

Application of two novel qPCR assays for quantifying pinniped environmental DNA (eDNA) in coastal environments

Julia Clem

PhD Candidate
University of Maine

3rd Marine Mammal eDNA Workshop
SMM 2024



Acknowledgements

Advisors and Co-Authors


- Kristina Cammen, University of Maine
- Michael Kinnison, University of Maine
- Lisa Sette, Center for Coastal Studies
- Christine Hudak, Center for Coastal Studies
- Justin Stevens, Maine Sea Grant
- Sydney Jackson, University of Maine
- Geneva York, University of Maine Environmental

CORE Facility

Funding

- Graduate stipend was provided by Maine-eDNA (NSF #1849227) and Maine Sea Grant
- Joshua A. Nickerson Conservation Fellowship
- Maine Outdoor Heritage Fund






Gray seal
Halichoerus grypus



Harbor seal
Phoca vitulina



Primer Development

in silico Design and Testing

- Identify gene region candidates
- Minimize potential for cross-amplification of non-target species

Specificity Testing

- Tissue extract from target and non-target species
- eDNA extract from rehabilitation pools

Assay Optimization

- Primer/probe concentrations
- Annealing temperature
- Amplification efficiency

Gray Seal Assay

106 bp
ND4 gene

Harbor Seal Assay

173 bp
cytb gene

Cape Cod, Massachusetts




Photo taken by the Center for Coastal Studies under NOAA Research Permit 21719.


Primer Deployment

Large gray seal aggregations (18-325 individuals)

Targeted sampling around haul-outs (along shoreline and 50m offshore)

Tested for gray seal

Penobscot River, Maine



Smaller seal aggregations (1-30 individuals)

Samples collected along a transect up the lower estuary

Tested for gray and harbor seals

Photo Credit: Lauri Leach

NOAA LOA 2016-22582

Cape Cod, Massachusetts




Photo taken by the Center for Coastal Studies under NOAA Research Permit 21719.



Detection Success Across Systems

92% samples collected 0 - 50 meters from shore detected gray seal eDNA at haulouts with 18-325 seals present



In 2 expanded surveys, gray seal eDNA was detected up to 150 meters from shore

Quantification in Space


Cape Cod,
Massachusetts



Quantification in Space



eDNA verses Abundance



Wide range of haul-out sizes offered the opportunity to test the influence of abundance on eDNA concentration

Together, distance from shore and seal abundance explained 73.88% of observed deviance in eDNA concentration

Distance had a stronger effect on eDNA concentration than seal abundance

Cape Cod, Massachusetts




Photo taken by the Center for Coastal Studies under NOAA Research Permit 21719.

Penobscot River, Maine




Photo Credit: Lauri Leach

NOAA LOA 2016-22582

Detection Success Across Systems

92% samples collected 0 - 50 meters from shore detected gray seal eDNA at haulouts with 18-325 seals present

In 2 expanded surveys, gray seal eDNA was detected up to 150 meters from shore

System characteristics and sampling scheme support a high rate of detection

Cape Cod, Massachusetts




Photo taken by the Center for Coastal Studies under NOAA Research Permit 21719.

Penobscot River, Maine




Photo Credit: Lauri Leach

NOAA LOA 2016-22582

Detection Success Across Systems

92% samples collected 0 - 50 meters from shore detected gray seal eDNA at haulouts with 18-325 seals present

In 2 expanded surveys, gray seal eDNA was detected up to 150 meters from shore

System characteristics and sampling scheme support a high rate of detection

Targeted analysis of samples collected during surveys with the highest number of seals observed (<30 individuals)

12% samples detected gray seal, 17% detected harbor seal

Cape Cod, Massachusetts




Photo taken by the Center for Coastal Studies under NOAA Research Permit 21719.

Penobscot River, Maine




Photo Credit: Lauri Leach

NOAA LOA 2016-22582

Detection Success Across Systems

92% samples collected 0 - 50 meters from shore detected gray seal eDNA at haulouts with 18-325 seals present

In 2 expanded surveys, gray seal eDNA was detected up to 150 meters from shore

System characteristics and sampling scheme support a high rate of detection

Targeted analysis of samples collected during surveys with the highest number of seals observed (<30 individuals)

12% samples detected gray seal, 17% detected harbor seal

Detection more sporadic, high rate of false positives


qPCR as a supporting data stream

Cape Cod, Massachusetts

- Accompanied by eDNA haplotyping
- qPCR data established a baseline understanding of eDNA signal in space around haul-out sites
 - Guide haplotyping efforts and help interpret detections

qPCR as a supporting data stream

Learn more about our eDNA
haplotyping efforts on
Monday!




Program Track 1

Ecology and Evolution - eDNA Application

 Mon, November 11

 11:30 AM - 12:30 PM

 [BelleVue Ballroom 2](#)

 Long and Short (or
Video) Presentations
Session

12:18 PM

Evaluating the Use of Environmental DNA for Pinniped Detection and Population Genetics at Haul-Out Sites

12:30 PM

Julia Sunnarborg¹, Christine Hudak², Lisa Sette², Geneva York³, Sydney Jackson⁴, Michael Kinnison⁴, Kristina Cammen⁴

¹University of Maine, Bangor, Maine, ²Center for Coastal Studies, Provincetown, Massachusetts,

³University of Maine Environmental DNA CORE Laboratory, Orono, Maine, ⁴University of Maine, Orono, Maine

qPCR as a supporting data stream


Cape Cod, Massachusetts

- Accompanied by eDNA haplotyping
- qPCR data established a baseline understanding of eDNA signal in space around haul-out sites
 - Guide haplotyping efforts and help interpret detections

Penobscot River, Maine




- Accompanied by 12S rRNA metabarcoding to characterize prey community

Metabarcoding – Where are the seals?



***No seal reads even
when visually observed***

Metabarcoding – Where are the seals?




Takeaways



qPCR assays created and field tested for gray and harbor seals



Conditions of the system and sampling setup influence detection rates and the potential for false negatives



Quantifying eDNA variation in space helps us understand the underlying processes of transport and signal dissipation for a given system



qPCR tools may provide a way to “pick out” seal detections where they might otherwise be swamped out

Thank you!

A group of seals, likely harbor seals, are resting in the shallow water. They are dark-colored with some lighter patches. Some seals are looking towards the camera, while others are resting their heads on the water. In the background, a small boat is visible on the left, and a sandy beach with more seals is in the distance under a clear blue sky.

Please reach out!
julia.sunnarborg@maine.edu